


SBRT Immobilization System

Meicen SBRT Immobilization System is highly accurate, comfortable, multi-functional and user-friendly device for SBRT treatment, which is compatible with the couchtops of Elekta, Varian or other brands. The paddle design on the belly and knee bridges can be adjusted in the direction of left and right, up and down to control the breathing range of the patient effectively during the treatment. The quick-lock and unlock system of the bridges is effortless. And, there are scales to record the fixed positions.

The system features a single-size breathing bridge that is fully adjustable to suit all patients. And the abundant scale marks improve the accuracy of positioning repeatability.

It can be used with Wing Board, Vacuum Cushions, Thermoplastic Masks, SRS Immobilization System and other products to provide various immobilization solutions for Radiation Therapy.



SBRT Immobilization System



High accuracy

Clear and precise scale, improves the accuracy of positioning repeatability.

Patient comfort


Ergonomically design for patient comfort.

Versatility

Multifunctional design, simple to form a set-up for various demand.

User-friendly

Easy and simple to operate.



SBRT Immobilization System


REF MC-STF001




Carbon Fiber SBRT Immobilization System Package

	Components	REF	Qty
MC-STF001 Standard Package	SBRT Platform	MC-STF001-A01	1
	Wing Board T-Grip	MC-WBF001	1
	Belly Bridge	MC-STF001-A02	1
	Knee Cushion	RE-FCA01	1
	Elevation Cushion	RE-FCA02	1
	Feet Cushion	RE-FCA03	1
	T-Shaped Vacuum Cushion	RC-E3010F	1
	SBRT Indexing Bar	MC-STF001-A07	4
Optional Product	Carbon Fiber Headrest	RE-HSF6	1
	Color Headrest	RE-HSB5	1
	Clear Headrest	RE-HSC6	1
	Knee Bridge	MC-STF001-A03	1
	SRS Cradle	MC-SRF001	1
	Indexing Bar (for Elekta couchtop)	RC-LBE02	2
	Indexing Bar (for Varian couchtop)	RC-LBV02	2
Shoulder Retractor	MC-STG001-A08	1	

Compatible with



S-type Violet Masks
page 36-39




Pin-lock type Violet Masks
page 46



SRS Immobilization System
page 8-11

SBRT Immobilization System

REF MC-STG001




MR SBRT Immobilization System Package

	Components	REF	Qty
MC-STG001 Standard Package	SBRT Platform	MC-STG001-A01	1
	Wing Board T-Grip	MC-WBG001	1
	Belly Bridge	MC-STG001-A02	1
	Knee Cushion	RE-FCA01	1
	Elevation Cushion	RE-FCA02	1
	Feet Cushion	RE-FCA03	1
	T-Shaped Vacuum Cushion	RM-E3010F	1
	SBRT Indexing Bars	MC-STG001-A07	4
Optional Product	Color Headrest	RE-HSB5	1
	Clear Headrest	RE-HSC6	1
	Knee Bridge	MC-STG001-A03	1
	SRS Cradle	MC-SRG001	1
	Indexing Bar (for Elekta couchtop)	RM-LBE01	2
	Indexing Bar (for Varian couchtop)	RM-LBV01	2
	Shoulder Retractor	MC-STG001-A08	1

SBRT Immobilization System - Components

SBRT Platform

MC-STF001-A01	Carbon Fiber
MC-STG001-A01	Glass Fiber MR



- Specification: 2000mm (L) x 530mm (W).
- Multiple corresponding holes, efficient indexing, install without restraint.
- Lightweight, easy to move.

Wing Board T-Grip (Page 28)

MC-WBF001	Carbon Fiber
MC-WBG001	Glass Fiber MR



- T-Grip design provides comfortable arms extension.
- Can be used with the T-Shaped Vacuum Cushion, RC-E3005F or RC-E3010F.

Belly Bridge

MC-STF001-A02	Carbon Fiber
MC-STG001-A02	Glass Fiber MR



- Ergonomic Compression Paddle for better comfort.
- The Paddle device can be adjusted in the direction of left and right, up and down, and has scales to record its fixed position.
- Abundant scale marks improve the accuracy of positioning repeatability.

Knee Bridge

MC-STF001-A03	Carbon Fiber
MC-STG001-A03	Glass Fiber MR




- The knee paddle device can be adjusted left and right, up and down, and has scales to record its fixed position.
- Abundant scale marks improve the accuracy of positioning repeatability.
- Knee paddle device is optional in double leg or single leg immobilization.

SBRT Immobilization System - Components

Knee & Feet Cushion (Page 61)

RE-FCA01	Knee Cushion
RE-FCA03	Feet Cushion
RE-FCA02	Elevation Cushion



Feet Cushion

Knee & Elevation Cushion

- Easy to operate, improve the comfort and stability of supine position during fixation.
- It can be directly fixed on the CT couchtop for use.

T-Shaped Vacuum Cushion (Page 54-56)

RC-E3005F	Non-open headrest
RC-E3010F	Open headrest



- Specification: 125x110cm.
- Open headrest area can be used together with the headrest to reduce the ray blocking.
- Provides arms-up positioning and support the torso region.


Shoulder Retractor MC-STG001-A08



- It can be used with the SBRT Platform to improve patient comfort effectively and reduce involuntary displacement.

Headrest (Page 58)

RE-HSB5	Color Headrest, 5 Sizes, A to E
RE-HSC6	Clear Headrest, 6 Sizes, A to F
RE-HSF6	Carbon Fiber Headrest, 6 Sizes, A to F




RE-HSB5

RE-HSC6

RE-HSF6

SBRT Indexing Bar

MC-STF001-A07	Carbon Fiber
MC-STG001-A07	Glass Fiber MR



- Allows positioning devices to be fixed to the SBRT Platform for enhancing reproducibility.